

Ephrin type-A receptor 4

Catalog No: #AP82969



Package Size: #AP82969-1 50ug #AP82969-2 100ug #AP82969-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ephrin type-A receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	EPH-like kinase 8,EK8,CEK8
Accession No.	Uniprot ID:Q07496
Uniprot	Q07496
GeneID	395559;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EPHA4

Note: This product is for in vitro research use only