

Poly(U)-specific endoribonuclease

Catalog No: #AP82971



Package Size: #AP82971-1 50ug #AP82971-2 100ug #AP82971-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Poly(U)-specific endoribonuclease
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Protein endoU,Uridylate-specific endoribonuclease
Accession No.	Uniprot ID:A6QLQ8
Uniprot	A6QLQ8
GeneID	525399;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ENDOU

Note: This product is for in vitro research use only