

Eukaryotic translation initiation factor 3 subunit A

Catalog No: #AP82973



Package Size: #AP82973-1 50ug #AP82973-2 100ug #AP82973-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Eukaryotic translation initiation factor 3 subunit A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	eIF3a, Centrosomin, Eukaryotic translation initiation factor 3 subunit 10, eIF-3-theta, eIF3 p167, eIF3 p180, eIF3 p185, p162, CsmA, Eif3, Eif3s10
Accession No.	Uniprot ID: P23116
Uniprot	P23116
GeneID	13669;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Eif3a

Note: This product is for in vitro research use only