

Ficolin-1

Catalog No: #AP82975



Package Size: #AP82975-1 50ug #AP82975-2 100ug #AP82975-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ficolin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Collagen/fibrinogen domain-containing protein 1,Ficolin-A,Ficolin-alpha,M-ficolin,FCNM
Accession No.	Uniprot ID:O00602
Uniprot	O00602
GeneID	2219;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FCN1

Note: This product is for in vitro research use only