

Dystonin

Catalog No: #AP82981



Package Size: #AP82981-1 50ug #AP82981-2 100ug #AP82981-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dystonin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Bullous pemphigoid antigen 1,BPA,Dystonia musculorum protein,Hemidesmosomal plaque protein,Microtubule actin cross-linking factor 2,Bpag1,Macf2
Accession No.	Uniprot ID:Q91ZU6
Uniprot	Q91ZU6
GeneID	13518;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dst

Note: This product is for in vitro research use only