

## Estrogen receptor beta

Catalog No: #AP82982



Package Size: #AP82982-1 50ug #AP82982-2 100ug #AP82982-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Estrogen receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ER-beta,Nuclear receptor subfamily 3 group A member 2,ESTRB,NR3A2
Accession No.	Uniprot ID:Q92731
Uniprot	Q92731
GeneID	2100;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ESR2

Note: This product is for in vitro research use only