

Sphingolipid delta(4)-desaturase DES1

Catalog No: #AP82984



Package Size: #AP82984-1 50ug #AP82984-2 100ug #AP82984-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sphingolipid delta(4)-desaturase DES1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Degenerative spermatocyte homolog 1, Degenerative spermatocyte-like protein RDES, Degs, Des1
Accession No.	Uniprot ID:Q5XIF5
Uniprot	Q5XIF5
GeneID	58970;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Degs1

Note: This product is for in vitro research use only