

Fatty acid desaturase 3

Catalog No: #AP82988



Package Size: #AP82988-1 50ug #AP82988-2 100ug #AP82988-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fatty acid desaturase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytochrome b5-related protein,CYB5RP
Accession No.	Uniprot ID:Q9Y5Q0
Uniprot	Q9Y5Q0
GeneID	3995;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FADS3

Note: This product is for in vitro research use only