

Protein FAM40A

Catalog No: #AP82992



Package Size: #AP82992-1 50ug #AP82992-2 100ug #AP82992-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein FAM40A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q0P5J8
Uniprot	Q0P5J8
GeneID	511120;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FAM40A

Note: This product is for in vitro research use only