

Growth arrest and DNA damage-inducible protein GADD45 gamma

Catalog No: #AP82996

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AP82996-1 50ug #AP82996-2 100ug #AP82996-3 1mg

Description

Product Name	Growth arrest and DNA damage-inducible protein GADD45 gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q9WTQ7
Uniprot	Q9WTQ7
GeneID	291005;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gadd45g

Note: This product is for in vitro research use only