

## Estrogen receptor

Catalog No: #AP82997



Package Size: #AP82997-1 50ug #AP82997-2 100ug #AP82997-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Estrogen receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ER,ER-alpha,Estradiol receptor,Nuclear receptor subfamily 3 group A member 1,Esr,Estr,Estra,Nr3a1
Accession No.	Uniprot ID:P19785
Uniprot	P19785
GeneID	13982;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Esr1

Note: This product is for in vitro research use only