

Desmoglein-1

Catalog No: #AP82999



Package Size: #AP82999-1 50ug #AP82999-2 100ug #AP82999-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Desmoglein-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Desmosomal glycoprotein 1,DG1,DGI
Accession No.	Uniprot ID:Q9GKQ8
Uniprot	Q9GKQ8
GeneID	403401;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DSG1

Note: This product is for in vitro research use only