

Desmoglein-2

Catalog No: #AP83000



Package Size: #AP83000-1 50ug #AP83000-2 100ug #AP83000-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Desmoglein-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cadherin family member 5,HDGC,CDHF5
Accession No.	Uniprot ID:Q14126
Uniprot	Q14126
GeneID	1829;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DSG2

Note: This product is for in vitro research use only