

# Protein FAM214A

Catalog No: #AP83009



Package Size: #AP83009-1 50ug #AP83009-2 100ug #AP83009-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein FAM214A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RCJMB04_30b20
Accession No.	Uniprot ID:Q5ZI58
Uniprot	Q5ZI58
GeneID	415416;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FAM214A

Note: This product is for in vitro research use only