

## G-protein coupled estrogen receptor 1

Catalog No: #AP83019



Package Size: #AP83019-1 50ug #AP83019-2 100ug #AP83019-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	G-protein coupled estrogen receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Chemoattractant receptor-like 2,G protein-coupled estrogen receptor 1,G-protein coupled receptor 30,G-protein coupled receptor 41,Membrane estrogen receptor,mER,Cmkr12,Gper,Gpr30,Gpr41
Accession No.	Uniprot ID:O08878
Uniprot	O08878
GeneID	171104;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gper1

Note: This product is for in vitro research use only