

## Dynein light chain 1, cytoplasmic

Catalog No: #AP83022

Package Size: #AP83022-1 50ug #AP83022-2 100ug #AP83022-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dynein light chain 1, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	8 kDa dynein light chain,DLC8,Dynein light chain LC8-type 1,Protein inhibitor of neuronal nitric oxide synthase,PIN,DLC1,DNCL1,DNCLC1,HDLC1
Accession No.	Uniprot ID:P63167
Uniprot	P63167
GeneID	8655;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:DYNLL1

Note: This product is for in vitro research use only