

## Estrogen-related receptor gamma

Catalog No: #AP83023



Package Size: #AP83023-1 50ug #AP83023-2 100ug #AP83023-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Estrogen-related receptor gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Estrogen receptor-related protein 3, Nuclear receptor subfamily 3 group B member 3, Err3, Kiaa0832, Nr3b3
Accession No.	Uniprot ID:P62509
Uniprot	P62509
GeneID	26381;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Esrrg

Note: This product is for in vitro research use only