

Gelsolin

Catalog No: #AP83026



Package Size: #AP83026-1 50ug #AP83026-2 100ug #AP83026-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gelsolin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Actin-depolymerizing factor,ADF,Brevin
Accession No.	Uniprot ID:Q68FP1
Uniprot	Q68FP1
GeneID	296654;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gsn

Note: This product is for in vitro research use only