

# High affinity immunoglobulin epsilon receptor subunit beta



Catalog No: #AP83031

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP83031-1 50ug #AP83031-2 100ug #AP83031-3 1mg

## Description

Product Name	High affinity immunoglobulin epsilon receptor subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FcERI,Fc epsilon receptor I beta-chain,IgE Fc receptor subunit beta,Membrane-spanning 4-domains subfamily A member 2,Fce1b,Fcer1b,Ms4a1
Accession No.	Uniprot ID:P20490
Uniprot	P20490
GeneID	14126;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ms4a2

Note: This product is for in vitro research use only