

Free fatty acid receptor 1

Catalog No: #AP83032



Package Size: #AP83032-1 50ug #AP83032-2 100ug #AP83032-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Free fatty acid receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	G-protein coupled receptor 40,Gpr40
Accession No.	Uniprot ID:Q8K3T4
Uniprot	Q8K3T4
GeneID	266607;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ffar1

Note: This product is for in vitro research use only