

Pikachurin

Catalog No: #AP83033



Package Size: #AP83033-1 50ug #AP83033-2 100ug #AP83033-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pikachurin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	EGF-like, fibronectin type-III and laminin G-like domain-containing protein,Nectican
Accession No.	Uniprot ID:Q4VBE4
Uniprot	Q4VBE4
GeneID	268780;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Egflam

Note: This product is for in vitro research use only