

GDP-L-fucose synthase

Catalog No: #AP83038



Package Size: #AP83038-1 50ug #AP83038-2 100ug #AP83038-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GDP-L-fucose synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GDP-4-keto-6-deoxy-D-mannose-3,5-epimerase-4-reductase,Protein FX,Red cell NADP(H)-binding protein,Short-chain dehydrogenase/reductase family 4E member 1,SDR4E1
Accession No.	Uniprot ID:Q13630
Uniprot	Q13630
GeneID	7264;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TSTA3

Note: This product is for in vitro research use only