

## GDP-L-fucose synthase

Catalog No: #AP83039



Package Size: #AP83039-1 50ug #AP83039-2 100ug #AP83039-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GDP-L-fucose synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GDP-4-keto-6-deoxy-D-mannose-3,5-epimerase-4-reductase,Protein FX,Red cell NADP(H)-binding protein,Transplantation antigen P35B,Tum-P35B antigen,P35b,Tstap35b
Accession No.	Uniprot ID:P23591
Uniprot	P23591
GeneID	22122;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tsta3

Note: This product is for in vitro research use only