

Fibroblast growth factor 22

Catalog No: #AP83040



Package Size: #AP83040-1 50ug #AP83040-2 100ug #AP83040-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 22
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FGF-22
Accession No.	Uniprot ID:Q9ESS2
Uniprot	Q9ESS2
GeneID	67112;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Fgf22

Note: This product is for in vitro research use only