

Fibroblast growth factor receptor 1

Catalog No: #AP83041

Package Size: #AP83041-1 50ug #AP83041-2 100ug #AP83041-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FGFR-1,Basic fibroblast growth factor receptor 1,BFGFR,Fms-like tyrosine kinase 2,FLT-2,N-sam,Proto-oncogene c-Fgr,CD331,BFGFR,CEK,FGFBR,FLG,FLT2,HBGFR
Accession No.	Uniprot ID:P11362
Uniprot	P11362
GeneID	2260;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGFR1

Note: This product is for in vitro research use only