

Fatty acid-binding protein, liver

Catalog No: #AP83045



Package Size: #AP83045-1 50ug #AP83045-2 100ug #AP83045-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fatty acid-binding protein, liver
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fatty acid-binding protein 1,Liver-type fatty acid-binding protein,L-FABP,FABPL
Accession No.	Uniprot ID:P07148
Uniprot	P07148
GeneID	2168;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FABP1

Note: This product is for in vitro research use only