

Fatty acid desaturase 2

Catalog No: #AP83047



Package Size: #AP83047-1 50ug #AP83047-2 100ug #AP83047-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fatty acid desaturase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Acyl-CoA 6-desaturase,Delta(6) fatty acid desaturase,D6D
Accession No.	Uniprot ID:O95864
Uniprot	O95864
GeneID	9415;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FADS2

Note: This product is for in vitro research use only