

Peptidyl-prolyl cis-trans isomerase FKBP8

Catalog No: #AP83055



Package Size: #AP83055-1 50ug #AP83055-2 100ug #AP83055-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peptidyl-prolyl cis-trans isomerase FKBP8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PPIase FKBP8,38 kDa FK506-binding protein,38 kDa FKBP,FK506-binding protein 8,FKBP-8,FKBPR38,Rotamase,Fkbp38,Sam11
Accession No.	Uniprot ID:O35465
Uniprot	O35465
GeneID	14232;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fkbp8

Note: This product is for in vitro research use only