

Coagulation factor V

Catalog No: #AP83060



Package Size: #AP83060-1 50ug #AP83060-2 100ug #AP83060-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Coagulation factor V
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Activated protein C cofactor
Accession No.	Uniprot ID:Q9GLP1
Uniprot	Q9GLP1
GeneID	397217;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: F5

Note: This product is for in vitro research use only