

Fibrinogen gamma-B chain

Catalog No: #AP83063



Package Size: #AP83063-1 50ug #AP83063-2 100ug #AP83063-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibrinogen gamma-B chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Gamma'
Accession No.	Uniprot ID:P12799
Uniprot	P12799
GeneID	280792;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FGG

Note: This product is for in vitro research use only