

## Fibroblast growth factor 2

Catalog No: #AP83066



Package Size: #AP83066-1 50ug #AP83066-2 100ug #AP83066-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fibroblast growth factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FGF-2,Basic fibroblast growth factor,bFGF,Heparin-binding growth factor 2,HBGF-2,FGFB
Accession No.	Uniprot ID:P09038
Uniprot	P09038
GeneID	2247;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FGF2

Note: This product is for in vitro research use only