

Fms-related tyrosine kinase 3 ligand

Catalog No: #AP83067



Package Size: #AP83067-1 50ug #AP83067-2 100ug #AP83067-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fms-related tyrosine kinase 3 ligand
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Flt3 ligand,SL cytokine
Accession No.	Uniprot ID:P49771
Uniprot	P49771
GeneID	2323;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FLT3LG

Note: This product is for in vitro research use only