

S-formylglutathione hydrolase

Catalog No: #AP83070



Package Size: #AP83070-1 50ug #AP83070-2 100ug #AP83070-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	S-formylglutathione hydrolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	FGH,Esterase D
Accession No.	Uniprot ID:Q9GJT2
Uniprot	Q9GJT2
GeneID	397127;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ESD

Note: This product is for in vitro research use only