

Ectonucleotide pyrophosphatase/phosphodiesterase family member 2

Catalog No: #AP83071

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83071-1 50ug #AP83071-2 100ug #AP83071-3 1mg

Description

Product Name	Ectonucleotide pyrophosphatase/phosphodiesterase family member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	E-NPP 2,Autotaxin,Extracellular lysophospholipase D,LysoPLD,Atx,Npps2
Accession No.	Uniprot ID:Q64610
Uniprot	Q64610
GeneID	84050;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Enpp2

Note: This product is for in vitro research use only