

## Gamma-enolase

Catalog No: #AP83072



Package Size: #AP83072-1 50ug #AP83072-2 100ug #AP83072-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gamma-enolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	2-phospho-D-glycerate hydro-lyase,Enolase 2,Neural enolase,Neuron-specific enolase,NSE,Eno-2
Accession No.	Uniprot ID:P07323
Uniprot	P07323
GeneID	100911625;24334;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Eno2

Note: This product is for in vitro research use only