

Epoxide hydrolase 4

Catalog No: #AP83073



Package Size: #AP83073-1 50ug #AP83073-2 100ug #AP83073-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Epoxide hydrolase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Abhydrolase domain-containing protein 7, Epoxide hydrolase-related protein, ABHD7, EPHXRP
Accession No.	Uniprot ID:Q8IUS5
Uniprot	Q8IUS5
GeneID	253152;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: EPHX4

Note: This product is for in vitro research use only