

## Alpha-2-HS-glycoprotein

Catalog No: #AP83074



Package Size: #AP83074-1 50ug #AP83074-2 100ug #AP83074-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Alpha-2-HS-glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	59 kDa bone sialic acid-containing protein,BSP,Fetuin-A,Glycoprotein PP63,Fetua
Accession No.	Uniprot ID:P24090
Uniprot	P24090
GeneID	25373;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ahsg

Note: This product is for in vitro research use only