

Fibroblast growth factor 1

Catalog No: #AP83075



Package Size: #AP83075-1 50ug #AP83075-2 100ug #AP83075-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	FGF-1,Acidic eye-derived growth factor II,EDGF II,Acidic fibroblast growth factor,aFGF,Endothelial cell growth factor,ECGF,Heparin-binding growth factor 1,HBGF-1,Prostatropin,AFGF,FGFA,HBGF-1
Accession No.	Uniprot ID:P03968
Uniprot	P03968
GeneID	281160;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGF1

Note: This product is for in vitro research use only