

# Fibroblast growth factor 5

Catalog No: #AP83077



Package Size: #AP83077-1 50ug #AP83077-2 100ug #AP83077-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fibroblast growth factor 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FGF-5,Heparin-binding growth factor 5,HBGF-5,Fgf-5
Accession No.	Uniprot ID:P15656
Uniprot	P15656
GeneID	14176;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Fgf5

Note: This product is for in vitro research use only