

Protein FAM203A

Catalog No: #AP83080



Package Size: #AP83080-1 50ug #AP83080-2 100ug #AP83080-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein FAM203A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Brain protein 16,Brp16,MNCb-5873
Accession No.	Uniprot ID:Q8C3I8
Uniprot	Q8C3I8
GeneID	59053;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fam203a

Note: This product is for in vitro research use only