

# Coagulation factor VIII

Catalog No: #AP83082



Package Size: #AP83082-1 50ug #AP83082-2 100ug #AP83082-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Coagulation factor VIII
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Procoagulant component,Cf8,F8c
Accession No.	Uniprot ID:Q06194
Uniprot	Q06194
GeneID	14069;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:F8

Note: This product is for in vitro research use only