

Lanosterol synthase

Catalog No: #AP83084



Package Size: #AP83084-1 50ug #AP83084-2 100ug #AP83084-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lanosterol synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	2,3-epoxysqualene--lanosterol cyclase,Oxidosqualene--lanosterol cyclase,OSC
Accession No.	Uniprot ID:Q8BLN5
Uniprot	Q8BLN5
GeneID	16987;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lss

Note: This product is for in vitro research use only