

DNA polymerase alpha subunit B

Catalog No: #AP83086



Package Size: #AP83086-1 50ug #AP83086-2 100ug #AP83086-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA polymerase alpha subunit B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DNA polymerase alpha 70 kDa subunit
Accession No.	Uniprot ID:P33611
Uniprot	P33611
GeneID	18969;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pola2

Note: This product is for in vitro research use only