

Fumarylacetoacetase

Catalog No: #AP83088



Package Size: #AP83088-1 50ug #AP83088-2 100ug #AP83088-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fumarylacetoacetase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FAA,Beta-diketonase,Fumarylacetoacetate hydrolase
Accession No.	Uniprot ID:P16930
Uniprot	P16930
GeneID	2184;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FAH

Note: This product is for in vitro research use only