

Fibrinogen C domain-containing protein 1

Catalog No: #AP83089



Package Size: #AP83089-1 50ug #AP83089-2 100ug #AP83089-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibrinogen C domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:A2AV25
Uniprot	A2AV25
GeneID	98970;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fibcd1

Note: This product is for in vitro research use only