

Eukaryotic translation initiation factor 1

Catalog No: #AP83092



Package Size: #AP83092-1 50ug #AP83092-2 100ug #AP83092-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Eukaryotic translation initiation factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	eIF1,Protein translation factor SUI1 homolog,SUI1
Accession No.	Uniprot ID:Q5E938
Uniprot	Q5E938
GeneID	509764;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EIF1

Note: This product is for in vitro research use only