

Galanin receptor type 2

Catalog No: #AP83097



Package Size: #AP83097-1 50ug #AP83097-2 100ug #AP83097-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galanin receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GAL2-R,GalnR2
Accession No.	Uniprot ID:O88854
Uniprot	O88854
GeneID	14428;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Galr2

Note: This product is for in vitro research use only