

CSK (Phospho-Tyr304) Antibody

Catalog No: #12494

Package Size: #12494-1 50ul #12494-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CSK (Phospho-Tyr304) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CSK (Phospho-Tyr304) Antibody detects endogenous levels of CSK only when phosphorylated at Tyr304
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CSK (Phospho-Tyr304)
Target Name	CSK
Modification	Phospho
Other Names	CSK, C-src tyrosine kinase, C-Src kinase, Protein-tyrosine kinase CYL, Tyrosine-protein kinase CSK, C-SRC, Cytoplasmic tyrosine kinase
Accession No.	Swiss-Prot#: P41240NCBI Gene ID: 1445
Uniprot	P41240
GeneID	1445;
Target Species	human
Calculated MW	55kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of MCF-7 cell (Left) and MCF-7 cell treated with 100ng/mL LPS for 30min (Right), using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only