

# Guanine nucleotide-binding protein G(s) subunit alpha isoforms XLas



Catalog No: #AP83102

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP83102-1 50ug #AP83102-2 100ug #AP83102-3 1mg

## Description

Product Name	Guanine nucleotide-binding protein G(s) subunit alpha isoforms XLas
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alternative gene product encoded by XL-exon,GNAS1
Accession No.	Uniprot ID:P84996
Uniprot	P84996
GeneID	2778;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:GNAS

Note: This product is for in vitro research use only