

Fatty acid-binding protein, brain

Catalog No: #AP83103



Package Size: #AP83103-1 50ug #AP83103-2 100ug #AP83103-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fatty acid-binding protein, brain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Brain lipid-binding protein,BLBP,Brain-type fatty acid-binding protein,B-FABP,Fatty acid-binding protein 7,Mammary-derived growth inhibitor related,BLBP,FABPB,MRG
Accession No.	Uniprot ID:O15540
Uniprot	O15540
GeneID	2173;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FABP7

Note: This product is for in vitro research use only